# ACCLEAN™

## ULTRASONIC SCALER

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by  HENRY SCHEIN®

## INDICATION FOR USE

This magnetostrictive ultrasonic scaler is designed for use during dental hygiene procedures and for periodontal debridement of periodontal diseases.

## HEALTH AND SAFETY

### Caution

- The Acclean Ultrasonic should not be used for restorative dental procedures involving the condensation of amalgam.
- Persons with cardiac pacemakers, defibrillators and other active implanted medical devices are cautioned that some types of electronic equipment might interfere with the operation of the device. Dental professionals should contact the device manufacturer or the patient's physician for specific recommendations.
- Do not use in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
- Do not modify the Acclean Ultrasonic. Modifications will invalidate the warranty and may endanger the patient and/or operator.

## GENERAL OPERATION/PRECAUTIONS

### Caution

- For professional use only.
- Use personal protective wear during use.
- Sterilize handpiece sheath and inserts prior to first use.
- Never use with inserts that show excessive wear.
- Close manually shut-off valve on the dental office water supply every night before leaving the office.
- Continuous water flow must be maintained at all times during operation to prevent overheating.
- Keep away from heat and direct sunlight.
- If overheating or any abnormalities occur, stop using immediately and contact Henry Schein® ProRepair.
- If excessive wear of insert is observed, discard and replace insert immediately.
- Bent, damaged or reshaped inserts are susceptible to breakage during use, discard these inserts immediately and replace.
- To prevent contact with the ultrasonic insert tip, be sure to retract patient's lips, cheeks, and tongue.
- Excessive heat may damage the system's electronics. Place the system where air is free to circulate on all sides and beneath it.
- The system is portable, but must be handled with care when moving.
- Equipment flushing and dental water supply system maintenance are strongly recommended.
- Always ensure that the electrical connections on the handpiece cable and the handpiece are clean and dry before use.

## PARTS AND CONTROLS

### IN THE BOX

1. Main Unit
2. Handpiece assembly
3. External power supply
4. Water supply tube with filter & 2 quick disconnects
5. Foot control pedal
6. Handpiece sheath (spare)

### MAIN UNIT CONTROLS

1. Power setting dial
2. Water flow adjustment
3. Boost indicator light
4. Cruise Mode Button
5. Purge Mode Button

### BACK OF UNIT

1. Power on/off switch
2. Foot pedal connection
3. Power supply connection
4. Water supply line connection

## SETUP

1. Place the Acclean Ultrasonic main unit on a level surface.
2. Confirm the power switch on the back of the main unit is in the OFF position.
3. Completely insert water supply tube into water supply line connection on the main unit (Fig. 1) until it snaps into place and then connect the quick disconnect (Fig. 2) to the water source. Confirm there are no leaks.

Fig. 1

**Note:** Only connect to a clean, filtered water supply capable of supplying 25 - 40 psi that has a manual shut-off valve so that water can be completely shut off when the office is unoccupied.

4. Firmly plug the foot control pedal connector to the foot pedal connection on the back of the unit.


Fig. 2

5. Connect handpiece assembly at the front of the unit and place handpiece in holder

**Note:** The red indicator on the handpiece assembly should not be seen if connected properly. Water leakage may occur if the handpiece assembly isn't properly connected to the unit. (Fig. 3)

Fig. 3

6. Connect the power cord from the external power supply into the back of the main unit then plug into power outlet.

# OPERATION

## ⚠ Caution


- Use personal protective wear during use.
- Do not use ultrasonic scaler if water is not continuously flowing to the insert tip.

## DAILY PREPARATION

1. Open manual shut-off valve on office water supply.  
**NOTE:** Remember to close valve at the end of the day.
2. Turn the system ON using the ON/OFF power switch at back of unit. Verify that Cruise and Purge buttons illuminate.
3. Set power setting dial to minimum setting.
4. Turn water flow adjustment dial to maximum setting.
5. Hold handpiece (without an insert installed) over sink or drain. Press the Purge button to activate Purge mode. The light will change from White to Green. Water line will be flushed automatically for 2 minutes. Once the cycle is completed, the light will change back to White. To deactivate Purge mode, press the Purge button again or press down on the foot pedal.
6. Set the water control knob as preferred for use on patient.
7. When ready, place a sterilized insert into the handpiece.

## SCALING

1. Slide sterilized autoclavable sheath onto handpiece and screw on clockwise until secure. Do not twist handpiece tubing and do not over tighten (Fig 4).

2. Wet o-ring on sterilized insert and install completely by twisting while pushing insert into handpiece until it snaps into place. Unit will automatically detect insert type.
3. Confirm power and water flow dials are to desired settings. Hold handpiece over sink and activate system by depressing foot control halfway. Verify water is flowing to end of insert tip.
4. Adjust water flow to ensure adequate cooling and irrigation for selected power setting.
5. To activate Boost feature depress foot control pedal completely; Boost indicator light will change from White to Orange during Boost mode. It may be necessary to adjust water flow during Boost for adequate cooling.
6. Cruise mode eliminates the need for the user to continually press on the foot pedal and can be activated by pressing the Cruise button on the unit. The light will change from White to Blue. Once activated, simply tap once on the foot pedal to begin scaling. To stop scaling, tap on the foot pedal again. Pressing the Cruise button again will de-activate Cruise mode and the light will revert back to White.
7. It is recommended to begin procedure using minimum power and to increase as necessary. Use a light touch while scaling and move handpiece in a circular or wave pattern.

### Caution

- Sterilize ultrasonic insert and handpiece sheath by autoclave before first use and after each patient.
- Do not immerse handpiece or sheath in liquid disinfectant.
- Always turn the device off via the switch and disconnect it from the electrical outlet before performing cleaning and disinfection.
- The device is not protected against penetration of liquids. Do not spray liquids directly onto the surface of the device.

### CLEANING & STERILIZATION

1. After each procedure, remove insert and sterilize according to manufacturer instructions.
2. Set water flow adjustment dial to maximum setting, hold handpiece over sink or drain, and activate Purge mode to flush water line.
3. Remove handpiece sheath (Fig 3).
4. Clean work area surfaces and Acclean Ultrasonic (including controls, cables, foot pedal, and tubing) with a non-immersion disinfecting solution. When cleaning the Acclean Ultrasonic, moisten towel with solution and wipe down unit. Discard towel in biohazard bin.  
**CAUTION:** Do not apply disinfecting solution directly to Acclean Ultrasonic.
5. Place handpiece sheath in a sterilization pouch and seal.
6. Steam autoclave for 20 minutes at 275° F (135° C).
7. Allow a minimum of 20 minutes for the drying cycle.

### WATER LINE FILTER MAINTENANCE

**NOTE:** It is recommended that you replace water filter every 3 months or when filter becomes discolored. Though your system comes with a water filter, it is recommended that you also use a filter in your dental water system supply line.

1. Confirm power switch on back of main unit is in the OFF position.
2. Close manual shut-off valve on office water supply.
3. Disconnect water supply tube from water source.
4. Grasp sleeve clamp on both sides of filter and twist counterclockwise to disconnect filter (Fig. 5).
5. Remove filter from both sides of water supply tube.
6. Install replacement filter onto water supply tube as shown in Figure 6.
7. Grasp sleeve clamp on both sides of filter and twist clockwise to hand tighten (Fig 4). Do not over tighten.
8. Reconnect water supply tube and follow instructions in Daily Preparation section to bleed air and test for leaks.


Fig. 5


Fig. 6

## ACCESSORIES

Acclean Ultrasonic Scaler

570-2773 Handpiece Assembly

570-2772 Handpiece Sheath

570-2774 Foot Pedal

570-0396 Replacement Water Filters (4 pack)

570-0397 Water Supply Tube with Filter and Quick Disconnect

Acclean Ultrasonic Inserts  
 112-6721 PRO IF-50 25k  
 112-6722 PRO-IF 50 30k  
 112-6723 PRO IF-1000 25k  
 112-6724 PRO IF-1000 30k  
 112-6725 ERGO IF-50 25k

112-6726 ERGO IF-50 30k  
 112-6727 ERGO IF-1000 25k  
 112-6728 ERGO IF-1000 30k  
 112-6733 X-Thin PRO 30k  
 112-6734 X-Thin ERGO 30k

112-6729 Classic P10 25k  
 112-6730 Classic P10 30k  
 112-6731 Classic P100 25k  
 112-6731 Classic P100 30k

## TROUBLESHOOTING

### AUDIBLE ALERTS

User Operation	Alert Description	Sound
Power Switch	Power ON	Beep x2
Cruise Button	Cruise Mode ON	Beep x1
Cruise Button	Cruise Mode OFF	Beep x2
Purge Button	Purge Mode ON	Beep x1
N/A	Purge cycle finished	Beep x2
Purge Button	Purge Mode OFF	Beep x2
Foot Pedal	Handpiece without insert	Beep x2, Pause, Beep x3
Foot Pedal	Foot Pedal does not connect with main unit	Beep x3

PROBLEM	CAUSE	SOLUTION
Unit will not operate; Touch Switch White LED light not on.	Unit not powered on.	Check AC power cord connections at Power Supply outlet Check unit DC power cable connections at back of unit. Make sure power outlet is live.
Unit will not operate; Touch Switch White LED light on.	Foot Pedal not securely connected.	Check the Foot Pedal connection at back of unit to confirm it's fully seated.
	Foot Pedal is damaged or defective.	Contact Henry Schein ProRepair.
Unit operates; but no water flow. Make sure water shut-off valve is open.	Poor water flow to/through unit.	Ensure water flow knob is properly adjusted.
	Insert is damaged or defective.	Replace insert.
	Water line filter is clogged.	Replace water line filter.
Unit operates; insert not vibrating.	Handpiece wire not securely connected.	Contact Henry Schein ProRepair.
	Insert not fully installed.	Install insert completely by twisting while pushing into handpiece until it snaps into place.

If the problem still can't be solved, please contact Henry Schein ProRepair.

## SPECIFICATIONS

Input Power: AC 100~240V~ 50/60Hz

Current: 1.4 A Max  
 Power: 60 W Max  
 Dimensions: 5.9" (L) x 4.9" (W) x 3.3" (H)  
 Weight: 1.3 lbs (Scaling Unit), 1.7 lbs (External Power Supply)  
 Tip Frequency: 25/30 kHz  
 Water Pressure: 25 - 40 psi

## SERVICE

For items that require service, including warranty service, please send to Henry Schein® ProRepair.  
 Henry Schein ProRepair  
 620 S. Placentia Ave.  
 Placentia, CA 92870 USA  
 Phone # 714 579-0175  
 Fax # 714 579-0186

## WARRANTY

This product is designed for use in a dental office and this warranty is not applicable to other uses. This product is warranted against manufacturing defects in materials and workmanship for 1 year from date of purchase.

Henry Schein® will repair or replace product at its sole discretion. This warranty shall be limited to replacement or repair of the product or its parts and shall not extend to any other claims, including but not limited to loss of profit, cost of removal or replacement, incidental or consequential damages or other similar claims arising from use of this product. Damage to product resulting from acts of God, faulty installation, misuse, tampering, accident, abuse, negligence, alterations, unauthorized repairs or problems unrelated to materials and workmanship are not covered by this warranty.

## INFORMATION ON ELECTROMAGNETIC COMPATIBILITY

Guidance and manufacturer's declaration-electromagnetic emissions		
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The product is suitable for use in locations in residential environments and in establishments directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuation / flicker emissions IEC 61000-3-3	Compliance	

## Guidance and manufacturer's declaration-electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	+ 2kV for power supply lines	+ 2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 1kV line(s) to line(s) + 2kV line(s) to earth	+ 1kV line(s) to line(s) + 2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT for 0,5 cycle 0% UT for 1 cycle 70% UT (30% dip in UT) for 25 cycles 0% UT for 5 s	0% UT for 0,5 cycle 0% UT for 1 cycle 70% UT (30% dip in UT) for 25 cycles 0% UT for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the product requires continued operation during power mains interruptions, it is recommended that the product be powered from an uninterruptible power supply or a battery.
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	The product power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**NOTE:** UT is the a.c. mains voltage prior to application of the test level.

### Guidance and manufacturer's declaration-electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below.  
The customer or the user of the product should assure that is used in such and environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHz to 80 MHz	3 Vrms 150 KHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the product including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. 
Radiated RF IEC 61000-4-3	6 V in ISM bands between 0,15 MHz and 80 MHz	6 V in ISM bands between 0,15 MHz and 80 MHz	
	80 % AM at 1 kHz  10 V/m 80MHz to 2,7 GHz	80 % AM at 1 kHz  10 V/m 80MHz to 2,7 GHz	

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the product is used exceeds the applicable RF compliance level above, the product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the product.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

### Recommended separation distance between portable and mobile RF communications equipment and the product

The product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,7 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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**Henry Schein Inc.**