# Maxima®

(EN) Electric System

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## INSTRUCTIONS FOR USE

Read before use. Keep in a safe place for future reference.

**⚠ CAUTION:** Federal (US) law restricts this device to sale by or on the order of a licensed professional.

### INDICATION FOR USE

The Maxima Electric System is intended to convert pneumatic output from a dental treatment center into electrical energy in order to drive an electric micromotor that operates electrical-ly-driven dental handpieces. This system is designed for use by a trained professional in the field of general dentistry.

### ELECTRICAL SAFETY INFORMATION

- Do not handle the power cord with wet hands to prevent an electric shock.
- Do not expose the device to water to prevent an electric shock.
- Do not use the device in a place where flammable materials are found due to risk of explosion.
- No part of the control unit shall be repaired by the user. Do not disassemble the device.
- Place the device on a stable location. Do not drop the device or exert an external force that may damage the device.
- Do not use a mobile phone or wireless RF communication device close by as such adjacent use may adversely influence this electrical device.

#### **⚠ Caution**

- This product shall be used only under the rated power specified. Power outside the rate indicated may damage unit
- Do not use device with a handpiece that does not meet the standard requirements.
- If motor begins to run rough, discontinue use and send in for request to the supplier.
- Check if the speed is set within the acceptable speed range prior to using the speed control button on the control panel.
- Use only the adapter that is provided by the manufacturer for this device.
- Do not disassemble or tamper with micromotor or control unit as it may damage the device.
- A Micromotor generates more power and torque than a traditional air turbine. Therefore, a handpiece that is equipped with micromotor that is poorly serviced or damaged can over-heat and cause serious burn injuries to patients.
- The following guidelines should be followed to ensure the safe operation of electric handpieces:
  - Follow maintenance instruction of handpiece manufacturer.
  - Examine the handpiece for damage before each use
  - Never apply pressure on the chuck button while handpiece is in operation or use hand-piece as a cheek or tongue retractor.

## GENERAL OPERATION / PRECAUTIONS

1. Read the instructions for use thoroughly prior to the first use.
2. Use the device only in accordance with the instructions for use.
3. Pay a close attention to safety of the patient while using the device.
4. Dispose the device and its parts in accordance with the applicable local laws and regulations.
5. For any issues relating to the use of the product including technical inquiries, please contact Henry Schein ProRepair.

## SPECIFICATIONS

Input	100-240V~, 1.4-0.7A, 47-63Hz
Power Consumption	105 VA(Max.)
Dimension	133 × 170 × 35 mm [Width × Length × Height]
Speed range	1,000 ~ 40,000 rpm
Gear ratio display	1:5, 1:1
Spray on/off	0
Program memory	3
Select a rotational direction	0
Class of protection	IPX0


## Micromotor

Max. Speed	40,000 rpm
Dimension	Ø20 × L63 mm
Weight	69g
Optic	White LED
Internal Irrigation	Yes
Coupling type	E-type connection per ISO 3964 type 3
Max. torque	3.0 N-cm
Air intake Max. quality per ISO 7494-2	4 bar (58 psi)
Water pressure intake Max.	2.5 bar (29 psi)

## AC/DC Adapter


Input	100-240 VAC , 1.4-0.7A
Frequency	47-63Hz
Output	36-38 VDC, 2.92 A
Power	105 W
Dimension	76 × 146 × 40 mm [Width × Length × Height]

## DISPLAY AND CONTROLS


NO.	DESCRIPTION	FUNCTION
1	Rotational Direction	Used to select a rotational direction (Fwd/Reverse) of the micromotor.
2	Rotational Direction Indicator	Displays the direction of the micromotor.
3	Micromotor LED ON/OFF Button	Used to turn on and off LED.
4	Micromotor LED Indicator	Indicates whether micromotor LED is running (LED On signals LED is turned on).
5	Spray ON/OFF Button	Used to switch on or off the spray function.
6	Spray Indicator	Displays the operational state of the spray function; water is sprayed when it is On.
7	Speed Controller	Used to set the operating speed of the motor.
8	Gear Ratio Selection Indicator	Displays the selected gear ratio of the handpiece.
9	Memory Program Selection Button	Used to save the program values set in the memory and select the saved program.
10	Gear Selection Button	Used to select the gear ratio of the handpiece connected to the motor.
11	Speed Indicator	Displays the operating speed of the handpiece.

## Rear

NO.	DESCRIPTION	FUNCTION
1	Tubing Connector-Standard 4-Hole Handpiece Hose Connector	Connects handpiece hose from dental unit to air & water supply tubing connector per ISO 9168 type 3 standard.
2	Motor Tubing / Cable	Supplies power and distributes air and water to micromotor.
3	Power Switch	Main power switch for unit.
4	Adapter Connector	Supplies power to unit.

## Micromotor

NO.	DESCRIPTION	FUNCTION
1	Handpiece coupling connector	Connect to attachment per ISO 3964 type 3 standard.
2	Body	Convert electrical energy to rotational motion.
3	LED	Provides light.

## Additional Device Information

1. Type of protections against electric shock: Class I equipment
2. Degree of protections against electric shock: Type B applied part
3. Mode of operation: Continuous operation
4. Applied parts  
Handpiece: Used in connection with a micromotor Handpiece sold separately.

## Environmental Conditions (Storage, Relocation, Operation, Shelf-life)

1. Storage condition  
Temperature: 0°C ~ +50°C      Humidity: 10 ~ 80%      Air pressure: 500hPa ~ 1060hPa
2. Relocation conditions  
Temperature: 0°C ~ +50°C      Humidity: 10 ~ 80%      Air pressure: 500hPa ~ 1060hPa
3. Operation conditions  
Temperature: +10°C ~ +35°C      Humidity: 30 ~ 80%      Air pressure: 700~1060 hPa
4. Shelf-life  
1 year based on standard fatigue shelf life of bearings the key components in the micromotor.

### ⚠ Caution


- Operation of the device in an environment other than those specified by manufacturer may damage the device.

## INSTALLATION & SET UP

### ⚠ Caution

- Check for damage to the cable and the connection pins prior to connecting with the micromotor.
- Check system connection diagram and cable / motor connection illustration below prior to installation.

# System Connection Diagram


1. Maxima Electric System
2. Maxima Electric Micromotor and hose
3. AC/DC Adapter 4-AC power plug

## Installing The Bracket

1. Hold metal bracket underneath the tray in the desired area and fasten with 6 appropriate sized screws.
2. Note: Using screws that are too big may damage bracket and tray. Using screws that are too small will not fasten the bracket to the tray.
3. Ensure a firm connection by lightly pressing down on the bracket with hand.
4. Slide Maxima Electric System onto the bracket until it settles into place.

## Connecting The Micromotor With The Cable / Tubing And Handpiece

1. Identify slot on micromotor and pin on motor tubing connector (Fig.1)
2. Align and insert the pin on motor tubing connector into slot on the motor terminal (Fig.2)
3. Turn the motor back nut clockwise to secure the connection. (Fig.3)
4. To attach a handpiece (handpiece sold separately) to the micromotor, identify guide stud on handpiece and hole on micromotor. (Fig. 4)
5. Lightly press the handpiece against the micromotor while turning it slightly counterclockwise until the guide stud locks into hole. Pull the handpiece back to make sure the handpiece is securely attached to the micromotor. (Fig.5) The micromotor coupling is designed to fit most E-Type electric handpieces for dental procedures per ISO 3964 type 3 standard.

## Handpiece Tubing Connector

1. Confirm that the delivery system tubing matches with the back terminal of the control unit for connection (Fig. 6) (coupling connection per ISO 9168 type 3 standard 4-hole handpiece hose)
2. Connect the tubing to the control unit and secure the connection to prevent the air or water leakage (Fig.7)

Fig. 6

Standard 4-hole handpiece tubing connector

Fig. 7


## Connecting The AC/DC Adapter Connection

Caution

- Only use power adapter supplied with this unit.

1. Insert the AC/DC adapter plug into the adapter connector (Fig. 8)

AC/DC adapter plug Fig. 8

## GENERAL OPERATION / PRECAUTIONS

### Speed Adjustment

1. The speed can be set up to a maximum of 40,000 rpm (1:5 : 200,000 rpm, 1:1 : 40,000rpm) by pressing the speed control button on the display.

### Setting A Rotational Direction

1. The direction of the micromotor can be set by pressing the rotational direction selection button on the control panel when the motor is not running.
2. When rotational direction indicator light is off, the motor is running in clockwise direction.

### Setting Optic LED On/Off

1. The LED installed in the micromotor can be turned on or off by pressing the OPTIC button on the control panel.

### Setting Indicator For Gear Ratios Of Handpiece

1. The gear ratio can be displayed as 1:5, 1:1 depending on the gear ratio of the handpiece selected with the use of the control panel.
2. The gear ratio can be adjusted to 1:5, 1:1 in a sequential order each time the operator presses GEAR button. The selected gear ratio is highlighted.

### Setting Water Spray

1. The water spray function can be turned on or off by pressing the SPRAY button on the display.

## Memory Program

1. How to save memory program : Customize settings of rotational direction, optic and so on;  
Press and hold one of the three program buttons P1, P2, or P3 for 3 seconds or longer.
2. How to retrieve saved program : Select the saved program by pressing P1, P2 or P3 button.

## CLEANING & MAINTENANCE

### Manual Cleaning

#### NOTE:

- The Micromotor is provided Non-Sterile
- Clean and disinfect the micromotor and tubing before first use and after each patient
- Sterilize the motor after manual cleaning and disinfection.

#### ⚠ Caution

- Routinely monitor dental unit water quality to ensure microbial counts are below the levels for drinking water recommended by your local agency. The U.S. Centers for Disease Control and Prevention (CDC) recommends a level of <500 CFU/ml of heterotrophic water bacteria.
- After each patient use, disinfect the tubing waterlines include the handpiece and micromotor waterlines by flushing with cleaning water to prevent cross-contamination according to handpiece waterlines disinfection procedure of main dental unit.

### Manual Internal Cleaning For Micromotor And Motor Tubing:

1. The Maxima Electric System internal water lines must be cleaned through the main dental unit by connecting the standard 4-hole handpiece connector and following the instructions for handpiece and hose waterlines disinfection of the main dental unit.  
(using the same procedure to disinfect waterlines of high speed handpieces on the main dental unit)
  - It is recommend to use Waterline Maintenance Tablets for main dental unit, it is the simplest way to maintain clean dental unit waterlines.
2. Do not remove attached handpiece or micromotor from the connecting cable during internal waterline cleaning and disinfection.
3. Hold the handpiece attached to micromotor over the cuspidor or over a watertight container on main dental unit with sufficient capacity and press the purging key on the main dental unit. The setting time to clean and disinfect waterlines is based on the main dental unit instructions.(Fig.9)


Fig. 9

### Manual external cleaning for micromotor and motor tubing:

Recommendation for external cleaning accessories:

- Clean cotton cloth
- Soft brush-Nylon Bristles

- Alcohol 70% Isopropyl
  - Protective glove uses during cleaning procedure.
1. For control box, micro motor cable / tubing (non removable) can be disinfected by wiping off using soft cotton cloth with isopropyl alcohol.
  2. For the micromotor, Disconnect the motor from the motor cable / tubing by turning motor back nut counter- clockwise and gently pull the motor out. (Fig.10)


Fig. 10


- Use air spray to purge the waterline & remove excessive moisture from the micromotor. (Fig. 11)
- Use soft brush with isopropyl alcohol to clean exterior surface then wipe off with soft cotton cloth.
- Visual inspection the cleanness of the micromotor prior to sterilization cycle.

Fig. 11

Fig. 12

## STERILIZATION

Caution

- Micromotor must be sterilized before first use and between each patient use by moist heat accordance with ISO 17665-1 in a steam sterilizer (Autoclave)

**Note:** After manual cleaning(internal & external cleaning), follow the instruction below to autoclave the micromotor.

### Autoclave using Air Removal Steam type Sterilizer (pre-vacuum)

1. Insert the motor into sterilization bag and load into autoclave. (We recommended using sterilization bag with at least one paper side.) The motor should be processed in the autoclave set for a minimum of 4 minutes at a sterilization temperature of 132°C / 270°F and a minimum dry time of 30 minutes.

### For sterilization of attachment handpieces:

1. Follow sterilization instructions provided by the manufacture of the handpiece.

**Note:** Replace the O-Ring if:

1. Water or air leaks occur.
2. The handpiece vibrates abnormally.
3. The supply of water or air cannot be stopped.
4. It is difficult to install or remove the handpiece.
5. There is no specific replacement period.

### Changing the O-Rings

To replace the motor O-Rings: (Fig.13)

## 4 O-RINGS


Fig. 13

1. Squeeze the O-Ring together between your thumb and forefinger to form a loop.
2. Pull the O-Ring off the motor.
3. Repeat steps 1 and 2 for the other three O-Rings.
4. Slide on the new O-Rings

**NOTE:** Always replace all four O-Rings at once to ensure a secure connection between the motor and attachment.

To replace the two small O-Rings at the motor connection back end: (Fig. 14)

1. Lift up the O-Rings with a pair of tweezers and remove them.
2. Carefully replace the two O-Rings with new ones.

Fig. 14

**NOTE:** The ISO 11498 standard stipulates a durability of at least 250 sterilization cycles. We recommend that you have the Maxima Electric Micromotor serviced after 500 sterilization cycles or one year, whichever comes first.

## TROUBLESHOOTING / ERROR MESSAGE REFERENCE GUIDE

ERROR CODE	STATUS	CAUSE OF ERROR	SOLUTION
E1	Micromotor connection error	Poor connection with cable / tubing	Visually inspect the micromotor connection.
E2	Micromotor defect	Poor micromotor connection or micromotor damage	Replace the micromotor if the micromotor connection shows no issues.
E3	Micromotor overload	Micromotor is being overloaded	Stop running the micromotor and restart after powered down for 3-5 minutes.

PROBLEM	CAUSE	SOLUTION
The micromotor is not running	Poor connection with AC/DC Adapter	Check the power connection.
	Poor micromotor connection	Check the micromotor connection.
	Control panel malfunction	Contact the manufacturer.
	Controller malfunction	Contact the manufacturer.
	Micromotor malfunction	Contact the manufacturer.
The micromotor speed can't get controlled.	Poor micromotor connection	Check the micromotor connection.
	Control panel malfunction	Contact the manufacturer.
Micromotor rotational direction cannot be changed.	Poor micromotor connection	Check the micromotor connection.
	Control panel malfunction	Contact the manufacturer.
Optic LED is not working.	Failure of LED inside the micromotor	Contact the manufacturer.

## ACCESSORIES

570-0685 Maxima Electric Micromotor

570-0686 AC/DC Adapter

570-1058 O-Ring Kit

## ELECTROMAGNETIC COMPATIBILITY

### Electromagnetic Emission

The product is suitable for use in a specific electromagnetic environment. The customer and/or the user of the product should assure that it is used in an electromagnetic environment as described below.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT GUIDANCE
RF-emission CISPR 11	Group 1	The product use RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF-emission CISPR 11	Class A	The product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purpose.
Harmonic emissions IEC 6100-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

## Electromagnetic Immunity

The product is suitable for use in a specific electromagnetic environment. The customer and/or the user of the product should assure that it is used in an electromagnetic environment as described below.

IMMUNITY TEST	EC 60601-LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT GUIDANCE
Electrostatic discharge(ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/bursts IEC61000-4-4	± 2kV for power supply lines ± 1kV for input/output lines	± 2kV for power supply lines ± 1kV for input/output lines	Mains power quality should be that of a typical commercial and/or hospital environment
Surge IEC61000-4-5	± 1kV differential mode ± 2kV common mode	± 1kV differential mode ± 2kV common mode	Mains power quality should be that of a typical commercial and/or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% DIP IN UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% DIP IN UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial and/or hospital environment. If the user of the product requires continued operation during power mains interruptions, it is recommended that the product be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

IMMUNITY TEST	EC 60601-LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT GUIDANCE
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz 3V/ m</p> <p>80 MHz to 2.5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: <math>d = 1.2\sqrt{P}</math></p> <p><math>d = 1.2\sqrt{P}</math> for 80 MHz to 800 MHz <math>d = 2.3\sqrt{P}</math> for 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in Watt (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level b in each frequency range</p> <p>Interference may occur in the vicinity of equipment marked with the symbol described lateral</p>

**Note 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, people and animals.

**a -** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered, if the measured field strength in the location in which the product is used exceeds the applicable RF compliance level above, the product should be observed, additional measures may be necessary, such as reorienting or relocating the product.

**b -** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended Separation Distances Between Portable And Mobile Hf- Communications Equipment And The Product

The product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer and/or the user of the product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product, according on output power and frequency of the communications equipment , as recommended in the following table.

Rated maximum output power of transmitter in watts (W)	Separation distance according to the frequency meters (m) of transmitter in meters (m)		
	150 kHz to 80 MHz d = $1.2\sqrt{P}$	80 MHz to 800 MHz d = $1.2\sqrt{P}$	800 MHz d = $2.3\sqrt{P}$ to 2.5 GHz
GHz			
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

## PRODUCT DISPOSAL GUIDELINE

### Disposal of Medical Devices

1. Follow your country specific laws, directives, standards and guidelines for the disposal of used electrical devices.
2. Ensure that the parts are not contaminated on disposal.

### Disposal Of The Packaging Material

1. All packaging materials have been selected considering environmental protection and can be recycled. Please send old/used packaging materials to the collection and reprocessing centers appropriate per local laws and regulations to minimize waste.

## WARRANTY

This product is designed for use in a dental office and this warranty is not applicable to other uses. This product is warranted against manufacturing defects in materials and workmanship for 1 year from date of purchase.

Henry Schein® will repair or replace product at its sole discretion. This warranty shall be limited to replacement or repair of the product or its parts and shall not extend to any other claims, including but not limited to loss of profit, cost of removal or replacement, incidental or consequential damages or other similar claims arising from use of this product. Damage to product resulting from acts of God, faulty installation, misuse, tampering, accident, abuse, negligence, alterations, unauthorized repairs or problems unrelated to materials and workmanship are not covered by this warranty.

Distributed by  
**HENRY SCHEIN INC.**  
 135 DURYEA ROAD  
 MELVILLE, NY 11747 USA

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