MAXIMA[™] 2000

(EN) High Output Curing Light

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INSTRUCTIONS FOR USE

Read before use. Keep in safe place for future reference.

 **Caution:** Federal (US) law restricts this device to sale by or on the order of a licensed professional.

Warning

- Heat development: As is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible damage. Therefore, this high-performance curing light must only be operated by trained professionals.

Caution

- Do not operate curing light without light guide in place.
- Never look directly into light source or direct it towards eyes.
- Wear protective eyewear and ensure glare shield and/or anti-glare cup is secured in place before use.
- Light guide can become hot during continuous use; avoid prolonged exposure to skin or oral soft tissue.
- Avoid direct contact of light guide with skin, oral tissue, and tooth surface.
- Before use, place orange anti-glare cup onto light guide to ensuring minimum distance of 2 mm between light guide and oral tissue or tooth surfaces.
- Do not directly touch teeth or oral surfaces with light guide.
- Always position light guide over the cure material when curing; do not direct light towards oral soft tissue.

Caution

- Sterilize light guide prior to first use.
- Use a disposable barrier sleeve over handpiece during use to prevent cross-contamination.
- To avoid damage to LED handpiece, do not drop or subject it to excessive shock.
- Keep away from water and other liquids.
- Do not use cold sterile or disinfectant on your handpiece or light guide.
- If charger base falls into water, unplug immediately. Do not reach into water or liquid.
- If either charger or LED handpiece come into contact with water or other liquids, stop using immediately and send to Henry Schein[®] ProRepair for evaluation.
- After continuous use, it is necessary to allow unit to cool for at least 30 seconds before beginning another cycle.
- Do not disassemble handpiece.
- Keep away from heat and direct sunlight.
- If any abnormalities are found, stop using immediately and contact Henry Schein[®] ProRepair.

INDICATION FOR USE

The Maxima® 2000 is a light emitting diode (LED) dental curing unit used for the polymerization of light-cured materials within the range of 360 ~ 490 nm. It is only to be used by trained dental professionals.

PARTS AND CONTROLS


LED Handpiece

1. 10 mm Black-coated Light Guide	
2. Nose Cone	
3. Display Screen	
4. MODE / ▲ Button	
5. SELECT / ▼ Button	
6. Serial Number	
7. Fan Vents	
8. Charging Contacts	
9. ⏻ / SET Trigger	

Charger Base and DC Adapter

1. Charger Base	
2. Charging Contacts	
3. Built-in Sensor	
4. Temperature Scale (°F / °C) Button	
5. Charging Indicator (orange)	
6. Ready Indicator (green)	
7. Display Panel	
8. DC Adapter	
9. DC Connection	
10. Serial Number	

Accessories

1. Glare Shield	
2. Anti-glare Cup	
3. Anti-glare Tube	


SETUP

1. Insert light guide into nose cone until it clicks securely into place. To remove, firmly pull light guide straight out (Fig. 1).
2. Slide glare shield over nose cone with thicker side of rubber ring towards handpiece until rubber ring seats completely in groove between nose cone and handpiece (Fig. 2).
3. Place orange anti-glare cup onto light guide and adjust according to preference (Fig. 3).
4. Connect DC adapter to DC connection on bottom of charger base and then plug into outlet (Fig. 4).

Note: When plugged in, display panel will brighten and show a reading of 0000 mW/cm² and the current temperature, and the green ready indicator will illuminate indicating base is ready to charge.

5. Place handpiece in charger base to begin charging. While charging, the orange charging indicator will illuminate and handpiece display screen will show a charging battery icon. The green ready indicator will illuminate when handpiece is fully charged.

Note: Allow handpiece to fully charge (minimum 4 hours) prior to first use.

OPERATION

Caution

- Wear protective eyewear and ensure glare shield and/or anti-glare cup is secured in place before use.
- Light guide can become hot during continuous use; avoid prolonged exposure to skin or oral soft tissue.
- Only use the 10 mm black-coated light guide (570-0254) approved for use with this product.
- Do not allow the device to come into direct contact with patients.

Curing Operation

Basic Operation

- If display panel is not already illuminated, press  / SET trigger once to turn on handpiece. An audible beep will sound and display screen will show curing mode and time from previous cycle.
- Press  / SET trigger to begin curing cycle. Two (2) beeps will countdown last 2 seconds of cycle and a 3rd beep will sound at end of cycle.
- Press  / SET trigger to stop curing during a cycle.

Curing Modes

- To toggle through curing mode, press MODE /  button until desired mode appears.
 - M1: Full Power Mode (1, 2, and 3 second options)
 - M2: Ramp-up Mode (5 and 10 second options)
 - M3: Pulse Mode (10, 20, and 30 second options)
- Press  / SET trigger to set desired mode.

Time Setting

- Press SELECT / ▼ to toggle through cycle time options.
- Press  / SET to begin curing cycle. Two (2) beeps will countdown last 2 seconds of cycle and a 3rd beep will sound at end of cycle.
- After a curing cycle is complete, pressing SELECT / ▼ will toggle through curing times for the last mode selected.
Note: Handpiece will automatically go to sleep after 30 seconds of inactivity. Press  / SET trigger to turn it back on.

Sensor Operation

Note: Charger base must be plugged in to use this feature.
Keep handpiece still while measuring as any movement will result in readout variations.

1. Press Temperature Scale button to toggle between °F and °C.
2. Select curing mode and time to be tested on handpiece (see Curing Operation).
3. Place light guide flat on built-in sensor.
4. Start curing cycle by pressing  / SET trigger.
5. Display panel will show light intensity and temperature.

Charging

Caution

- Ensure handpiece is completely dry before placing in charger.
- Charge and store handpiece in cool environment. Keep away from heat sources.

Note: To maximize the life of your high performance lithium battery:

- Avoid overcharging; do not store handpiece on plugged-in charger base for extended periods of time.
- Do not store curing light when battery is fully charged or when nearly empty.
- Use curing light at least once a month.
- Allow handpiece to fully charge (minimum 4 hours) prior to first use.

1. Connect DC adapter to DC connection on bottom of charger base (Fig. 4) and then plug into outlet.
2. Place handpiece in charger base to begin charging. While charging, orange charging indicator will illuminate and handpiece display screen will show a charging battery icon. The green ready indicator will illuminate when handpiece is fully charged.

Note: If battery capacity is considerably reduced or will not hold a full charge, the battery needs to be replaced. Contact Henry Schein® ProRepair for repair service.

Charger Base Indicator Lights

- Orange charging indicator: Charging
- Green ready indicator when handpiece in base: Fully charged
- Green ready indicator when no handpiece in base: Ready to charge

Displays

Handpiece Display Screen

M1  Full Power Mode

M2  Ramp-up Mode

M3  Pulse Mode

Full wavelength active (360 ~ 490 nm)

Fully Charged Battery

Low Battery: Recharge

Charging Battery

Overheating: Handpiece will power off and can restart after cooling down. This may occur during long periods of high use.

Note: Cooling System

- This device is designed with a built-in cooling fan which automatically activates whenever device is in use and will continue to run for a few seconds after cycle has ended.
- If internal temperature exceeds 113° F / 45° C, fan will continue to run until internal temperature has cooled down.
- If internal temperature exceeds 149° F / 65° C, the device will automatically shut off and fan will continue to run. When handpiece has cooled down, it can be turned on again.

Charger Base Display Panel

1. Light intensity
2. Measuring full wavelength intensity
3. Temperature of light guide at tip
(for reference only)

CLEANING AND STERILIZATION

Caution

- Sterilize light guide before first use.
- Sterilize light guide after each patient use.
- Keep light guide in sterilization pouch until ready to use.
- Use a disposable barrier sleeve over handpiece during use to prevent cross-contamination.
- Wear appropriate personal protective equipment during the cleaning and sterilization process. This should include eye protection and gloves.
- Do not run or submerge handpiece in water.
- Do not use cold sterile or disinfectant on your handpiece or light guide.
- Do not sterilize handpiece body.

1. Remove anti-glare cup from light guide then remove light guide and glare shield from handpiece.
2. For 1 minute wipe handpiece, light guide, anti-glare cup, and glare shield with a damp cloth moistened with a 70% isopropyl alcohol solution.
3. Wipe everything again for 1 minute with a soft cloth soaked in warm water.
4. Place light guide in sterilization pouch and seal. Bag light guides separately.
5. Sterilize light guide at 270° F / 132° C at 75 psi for 4 min.

NOTE: Ensure that all crevices are thoroughly cleaned.

ACCESSORIES

570-0254 10 mm Replacement Light Guide

TROUBLESHOOTING

Problem	Cause	Solution
Curing light doesn't turn on	Battery not charged	Charge battery
	Battery or PC board failure	Contact Henry Schein ProRepair
Curing light doesn't charge	Power adapter not fully plugged in	Reconnect power adapter to charger base and check connection to outlet
	Charging contacts on handpiece not properly touching contacts in charger base	Properly seat handpiece in charger base
	Charging contacts dirty	Clean charging contacts on handpiece and charger base with dry soft cloth
	Charger stand or power adapter failure	Contact Henry Schein ProRepair
	Handpiece battery failure	Contact Henry Schein ProRepair
Handpiece display screen doesn't turn on	Battery not charged	Charge
	Battery or PC board failure	Contact Henry Schein ProRepair
Charger base display panel doesn't turn on	Power adapter not fully plugged in	Reconnect power adapter to charger base and check connection to outlet
	Charger stand or power adapter failure	Contact Henry Schein ProRepair

SPECIFICATIONS

LED HANDPIECE	
LED lamp	5 W
Light intensity:	M1 (Full Power): 2,000 mW/cm ² (min)
	M2 (Ramp-up): 1,600 mW/cm ²
	M3 (Pulse): 1,600 mW/cm ² (1 sec - strong) / 600 mW/cm ² (0.5 sec - weak)
Wavelength:	360 ~ 490 nm
Weight:	215 g
Dimensions:	250 x 40 x 90 mm
Battery:	7.4 V, 1350 mAh
Light guide:	75 x Ø10 → Ø10 mm (included)
CHARGER BASE	
Input power:	100~240V ac 50/60 Hz
Output power:	DC 9 V, 1.3 A
Weight:	260 g
Dimensions:	150 x 90 x 120 mm

GENERAL

System operational conditions	Temp: 77° ~ 95° F (25° ~ 35° C)
	Humidity: 15% - 95% at 104° F / 40° C, noncondensing
	Pressure: 700 to 1060 hPa
System storage conditions	Temp: -4° ~ 149° F (-20° ~ 65° C)
	Humidity: 15% - 95% at 104° F / 40° C, noncondensing
	Pressure: 700 to 1060 hPa
System transport conditions	Temp: -4° ~ 149° F (-20° ~ 65° C)
	Humidity: 15% - 95% at 140° F / 40° C, noncondensing
	Pressure: 700 to 1060 hPa

SERVICE

For items that require service, including warranty service, please send to Henry Schein ProRepair.

Henry Schein ProRepair
620 S. Placentia Ave.
Placentia, CA 92870 USA
Phone # 714-579-0175
Fax # 714-579-0186

WARRANTY

This product is designed for use in a dental office and this warranty is not applicable to other uses. This product is warranted against manufacturing defects in materials and workmanship 1 year from date of purchase. Henry Schein® will repair or replace product at its sole discretion. This warranty shall be limited to replacement or repair of the product or its parts and shall not extend to any other claims, including but not limited to loss of profit, cost of removal or replacement, incidental or consequential damages or other similar claims arising from use of this product. Damage to product resulting from acts of God, faulty installation, misuse, tampering, accident, abuse, negligence, alterations, unauthorized repairs or problems unrelated to materials and workmanship are not covered by this warranty.

DISPOSAL OF MEDICAL DEVICES

Any medical devices which could possibly be contaminated must first be decontaminated by the responsible doctor or medical institution and then be disposed of by an agent licensed and qualified to handle medical and industrial waste.

INFORMATION ON ELECTROMAGNETIC COMPATIBILITY

Guidance and manufacturer's declaration – electromagnetic emissions		
The Maxima® 2000 is intended for use in the electromagnetic environment specified below. The customer or the user of the Maxima® 2000 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Maxima® 2000 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Maxima® 2000 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Recommended separation distances between portable and mobile RF communications equipment and the Maxima® 2000			
The Maxima® 2000 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Maxima® 2000 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Maxima® 2000 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guidance and manufacturer's declaration – electromagnetic immunity

The Maxima® 2000 is intended for use in the electromagnetic environment specified below. The customer or the user of the Maxima® 2000 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % <i>UT</i> (>95% dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95% dip in <i>UT</i>) for 5 sec	<5 % <i>UT</i> (>95% dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95% dip in <i>UT</i>) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Maxima® 2000 requires continued operation during power mains interruptions, it is recommended that the Maxima® 2000 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: *UT* is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The Maxima® 2000 is intended for use in the electromagnetic environment specified below.
The customer or the user of the Maxima® 2000 should assure that it is used in such an environment

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	3 Vrms 3 A/m	Portable and mobile RF communications equipment should be used no closer to any part of the Maxima® 2000, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2 \sqrt{P}$ ----- $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz ----- $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Maxima® 2000 is used exceeds the applicable RF compliance level above, the Maxima® 2000 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Maxima® 2000.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than V/m.

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