


FDA Cleared

enbio.



Enbio^S
World's fastest autoclave.
Vacuum+drying.

Designed  in Switzerland

Bone dry pouches in no time.



15 min



7 min

Enbio can sterilize a load of wrapped instruments in just 15 minutes.

Full drying means that pouches are bone dry. Not somewhat dry.

In case of emergency, when you need an instrument sterilized for immediate use, FAST process does the job in just 7 minutes.

4 full cycles with drying per hour.



Enbio S is 6 times smaller and 3 times lighter than a typical vacuum autoclave. But when it comes to work efficiency, it keeps up the pace.


	Pouches / h	Instruments / h
Enbio ^S	32	50



The new definition of sterilization speed.



Time is a relative value. It flows differently under different circumstances.
But when you need to have your instruments ready, you need them, well, FAST.
While other autoclaves might test your patience, Enbio bends the laws
of physics to cut the waiting to the invisible.



Inaudibly silent.



Enbio autoclaves work unbelievably silent.
Wherever you use it, nobody will notice it's working.

Enbio^S. Speedmaster.

Enbio S is compact to fit even a confined space. Unbelievably short process times will save your procedure in emergency, or allow sterilization of the instruments as you go, between patients.




Vacuum
Pouches
+drying

FDA Cleared

<p>Unwrapped FAST / 134°C</p> <p>7 min</p>	<p>Wrapped Vacuum / 134°C</p> <p>15 min</p>	<p>Wrapped Vacuum / 121°C</p> <p>45 min</p>
---	--	--



reddot design award
winner



Unit weight
33 lb

Capacity
Max. packages amount (circa) **8**
Max. recommended batch weight **1,1 lb**
Process chamber volume **2,7 l**

Pressure & temperature
Working pressure **2,1 Bar**
Max. pressure **2,3 Bar**
Max. process temperature **137°C**


Dimensions
Width **9.8 in**
Height **6.3 in**
Depth **22 in**

Electrical
Power supply **110V/60Hz**
Installed power **3,23 kW**
Max. power consumption **15 A**

Misc
Ingress protection level **IP20**
Data archiving **USB stick**
Water conductivity* **< 15 µS/cm**

* - demineralized water only


Enbio user's manual



1. Plug in
2. Choose program
3. Start

Enbio Care

Perfect technology meets perfect care.
All Enbio devices are covered with world's unique
warranty and care program, with backup device
delivered to you for the time of repair.



20.000 units in 4 years.

In just 4 years, Enbio's worldwide community has grown to 20,000 users, who sterilize instruments in our autoclaves, for the complete safety of their customers and patients.





Designed
in Switzerland

Worldwide presence.

Enbio autoclaves are manufactured in the EU and distributed globally through logistics centers covering Europe, North America and the Middle East.



Enbio Group AG
Eichengasse 3, 4702 Oensingen
Switzerland

FDA Establishment Registration
Number: 3016777357

US Distribution Center
Enbio Corp.
3535 Gravel Springs Rd, Ext. 205
Buford, GA 30519